

# LED4500

Light. Years ahead.



**SUNOPTIC**

# The Economy of LED. The Power of Xenon.

Our most advanced surgical light source yet provides many reliable hours of pure white light with 20% more output than before.

## Key Features

### Optimized for Use with

- Surgical headlamps
- Endoscopic systems
- Lighted Instruments

### System Highlights

- 105,000 lux at 16"/40cm working distance\*
- 6500° color temperature for accurate tissue rendition
- Stand-by mode for ease of use
- Intuitive user interface
- Four port turret for compatibility with ACMI, Wolf, Olympus, and Storz connectors
- Bright, evenly distributed light spot with headlamps

### Sturdy and Reliable

- No lamp to replace
- Consistent illumination throughout product life
- Low-profile and whisper quiet operation
- Minimal service requirements
- Backed by our robust three year warranty



Shown with our SG40-110-VHB Headlight

\* Industry standard distance to the surgical site

*Designed, Developed, and Manufactured in the United States.*

